

# LeiLOX TPLO Swing



## Titanium for Best Biocompatibility

LeiLOX TPLO Swing implants are made of medical grade Titanium, which makes them very durable but lightweight, highly biocompatible, and low in temperature sensitiveness.



## Multiaxial Locking

The screws can be locked in a 90° angle with a 12° deviation in any direction. This allows you to angle the screws away from vital structures. The robust LeiStar screw head can be locked firmly into the plate.



## Anatomically Shaped Limited Contact Dynamic Compression Plates

The LeiLOX TPLO Swing plates are contoured to match the anatomic shape of the bone. This makes it easier to place the plate in an optimal position.

The plate features limited contact dynamic compression to minimize vascular damage to the plated bone segment.



## Double Compression

Two precisely designed compression holes enable a very tight compression of the osteotomy.



## Interchangeable 2.7 & 3.5mm and 2.0 & 2.4mm Screws

Because the screwheads are identical, all Titanium 2.7/3.5 LeiLOX plates (TPLO Swing, CBLO, and CCWO systems) work with 2.7mm as well as 3.5mm screws in all of the plate sizes. Same applies for the 2.0/2.4 systems. This offers flexibility and ideal implant selection for each patient. Moreover, this saves on inventory cost.

High Performance Sawblades are available. The Titanium Nitrate coating allows them to last significantly longer than standard blades.

# LeiLOX

## TPLO Swing

TPLO Swing  
Titanium



Special Thanks to  
Dr. Matan Or

# LeiLOX TPLO SWING

**R** RITA  
LEIBINGER  
MEDICAL

**2.0 / 2.4 TPLO Swing Set**

Contains:


- 1 TPLO Implants and Instruments Tray with Lid
- 2 of each TPLO Swing Locking Plate
- 3 of each Cortical Screw
- (2.0, 8-18mm // 2.4, 8-22mm, 42 total)
- 5 of each Locking Screw
- (2.0, 6-24mm // 2.4, 6-30mm, 115 total)
- 2 Drills (1.5 & 1.8mm)
- 5 K-Wires
- 2x2 Locking Drill Guides
- 2 Compression Drill Guides
- 1 TPLO Jig
- 1 Screwdriver Handle
- 1 Screwdriver Shaft T8
- 1 Depth Gauge



**142-2820-24**

Tray without contents

**142-2800-10**



**Sterilization Container**

See Sterilization Containers for Implants and Instruments Trays on **Page 204**.

**2.0 / 2.4 LeiLOX TPLO Swing Locking Plates Titanium**



LeiLOX TPLO Swing Locking Plate  
2.0 mm, left, 33 mm, Titanium

**142-2820-10**



LeiLOX TPLO Swing Locking Plate  
2.0 mm, right, 33 mm, Titanium

**142-2820-00**



LeiLOX TPLO Swing Locking Plate  
2.4 mm, left, 37 mm, Titanium

**142-2824-10**



LeiLOX TPLO Swing Locking Plate  
2.4 mm, right, 37 mm, Titanium

**142-2824-00**

## TPLO Screws Locking & Non-Locking 2.0/2.4 Titanium

### 2.0 / 2.4 LeiLOX Locking Screw Titanium

For LeiLOX Locking Systems, LeiStar T8  
self-holding (T8 Shaft from Rita Leibinger recommended)  
self-tapping with three flute cutting edge



	2.0 mm	2.4 mm
Length (mm)	Product Code	Product Code
06	242-220-06	242-224-06
08	242-220-08	242-224-08
10	242-220-10	242-224-10
12	242-220-12	242-224-12
14	242-220-14	242-224-14
16	242-220-16	242-224-16
18	242-220-18	242-224-18
20	242-220-20	242-224-20
22	242-220-22	242-224-22
24	242-220-24	242-224-24
26	242-220-26	242-224-26
28	242-220-28	242-224-28
30	242-220-30	242-224-30
32		242-224-32
34		242-224-34
36		242-224-36
38		242-224-38
40		242-224-40

### 2.0 / 2.4 Cortical Screw (Non-Locking) Titanium

LeiStar T8, Non-Locking  
self-holding (T8 Shaft from Rita Leibinger recommended)  
self-tapping with three flute cutting edge



	2.0 mm	2.4 mm
Length (mm)	Product Code	Product Code
08	245-520-08	245-524-08
10	245-520-10	245-524-10
12	245-520-12	245-524-12
14	245-520-14	245-524-14
16	245-520-16	245-524-16
18	245-520-18	245-524-18
20		245-524-20
22		245-524-22

See more screw lengths on **Page 126.**

Screw Racks on **Page 198.**



**2.7 / 3.5 TPLO Swing Set**

Contains:

- 1 TPLO Implant Tray with Lid
- 2 of each TPLO Swing Plate
- 3 of each Cortical Screw (16-34mm, 60 total)
- 5 of each Locking Screw (16-46mm, 160 total)

**142-2827-35**

Tray without contents

**142-2800-20**



**2.7 / 3.5 LeiLOX TPLO Swing Locking Plates Titanium**



LeiLOX TPLO Swing Locking Plate  
2.7 mm, left, 44 mm, Titanium

**142-2827-10**



LeiLOX TPLO Swing Locking Plate  
2.7 mm, right, 44 mm, Titanium

**142-2827-00**



LeiLOX TPLO Swing Locking Plate  
3.5 mm, left, 63 mm, Titanium

**142-2835-10**



LeiLOX TPLO Swing Locking Plate  
3.5 mm, right, 63 mm, Titanium

**142-2835-00**



LeiLOX TPLO Swing Locking Plate  
3.5 mm, left, 74 mm, Titanium

**142-2835-60**



LeiLOX TPLO Swing Locking Plate  
3.5 mm, right, 74 mm, Titanium

**142-2835-50**

# TPLO Screws Locking & Non-Locking 2.7/3.5 Titanium

## 2.7 / 3.5 LeiLOX Locking Screw Titanium

For LeiLOX Locking Systems, LeiStar T10 self-holding (T10 Shaft from Rita Leibinger recommended) self-tapping with three flute cutting edge



	2.7 mm	3.5 mm
Length (mm)	Product Code	Product Code
10	242-227-10	242-235-10
12	242-227-12	242-235-12
14	242-227-14	242-235-14
16	242-227-16	242-235-16
18	242-227-18	242-235-18
20	242-227-20	242-235-20
22	242-227-22	242-235-22
24	242-227-24	242-235-24
26	242-227-26	242-235-26
28	242-227-28	242-235-28
30	242-227-30	242-235-30
32	242-227-32	242-235-32
34	242-227-34	242-235-34
36	242-227-36	242-235-36
38	242-227-38	242-235-38
40	242-227-40	242-235-40
42	242-227-42	242-235-42
44	242-227-44	242-235-44
46	242-227-46	242-235-46
48	242-227-48	242-235-48
50	242-227-50	242-235-50
52		242-235-52
54		242-235-54
56		242-235-56
58		242-235-58
60		242-235-60

See more screw lengths on **Page 126.**

## 2.7 / 3.5 Cortical Screw Titanium

LeiStar T10, Non-Locking self-holding (T10 Shaft from Rita Leibinger recommended) self-tapping with three flute cutting edge



	2.7 mm	3.5 mm
Length (mm)	Product Code	Product Code
16	245-527-16	245-535-16
18	245-527-18	245-535-18
20	245-527-20	245-535-20
22	245-527-22	245-535-22
24	245-527-24	245-535-24
26	245-527-26	245-535-26
28	245-527-28	245-535-28
30	245-527-30	245-535-30
32	245-527-32	245-535-32
34	245-527-34	245-535-34




### TPLO Saw Blades

See radial saw blades for TPLO and CBLO in various radii on **Page 191.**




### Screw Racks

See screw racks for 1.5mm up to 3.5mm screws on **Page 198.**



### TPLO Instruments

See essential instruments for a TPLO Surgery on **Page 78.**













### Instrument Set

See Instrument Set for TPLO / CBLO / CCWO on **Page 78.**

**TPLO Patella Luxation Spacers**

For LeiLOX TPLO Stainless Steel and TPLO Swing Titanium Systems

Description	Compatibility	Product Code
<p>2.0 TPLO Patella Luxation Spacer Titanium, 2.0mm thick for 2.0 mm TPLO Locking Plates, left &amp; right</p> 		<b>142-8020-00</b>
<p>2.4 TPLO Patella Luxation Spacer Titanium, 2.0mm thick for 2.4 mm TPLO Locking Plates, left &amp; right</p> 		<b>142-8024-00</b>
<p>2.7 TPLO Patella Luxation Spacer Titanium, 2.0mm thick for 2.7 mm TPLO Locking Plates, left &amp; right</p> 		<b>142-8027-00</b>
<p>3.5 TPLO Patella Luxation Spacer Titanium, 2.0mm thick for 3.5 mm TPLO Locking Plates, left &amp; right</p> 		<b>142-8035-00</b>
<p>3.5 broad TPLO Patella Luxation Spacer Titanium, 2.0mm thick for 3.5 mm broad TPLO Locking Plates, left &amp; right</p> 		<b>142-8035-10</b>

**Tibia Tappet**

For inserting TPLO Patella Luxation Spacers

Product Code	Description
<b>142-2890-00</b>	For 2.0 / 2.4 mm plates
<b>142-2891-00</b>	For 2.7/3.5 mm plates

