



EQUINE DENTAL

R RITA
LEIBINGER
MEDICAL

Equine Dental Plug Set

NEW

Titanium Equine Dental Plugs

Expertly designed to close non-healing oro-nasal and oro-sinusal fistulae in horses, these plugs provide convenient insertion through a straightforward transbuccal approach using a screwdriver. Their threaded design ensures secure and stable placement. Crafted from high-quality titanium, the RITA LEIBINGER Dental Plugs are highly biocompatible and durable.



Equine Dental Plug Set


Set contains:

Product Code	Description	Qty
888-9500-10	Tray with Lid	1
123-3019-00	Screwdriver Shaft T30	1
123-1019-00	Screwdriver Shaft T10	1
247-0101-00	Silicone T-Handle, AO	1
247-0102-00	Screwdriver Handle, AO	1
888-9501-00	Screws for Dental Plugs, Titanium	10
888-9507-20	Equine Dental Plug, Ø7 x 20mm	2
888-9510-20	Equine Dental Plug, Ø10 x 20mm	2
888-9513-20	Equine Dental Plug, Ø13 x 20mm	2
888-9516-20	Equine Dental Plug, Ø16 x 20mm	2
888-9519-20	Equine Dental Plug, Ø19 x 20mm	2

888-9500-00

Tray without contents

888-9500-10



Sterilization Container


See Sterilization Containers for Implants and Instruments Trays on **Page 204**.

Equine Dental Plugs

Screwdriver Shafts for Dental Plugs

Equine Dental Plugs

Titanium, LeiStar, including screws



Product Code	Ø	Length
888-9507-20	7 mm	20 mm
888-9510-20	10 mm	20 mm
888-9513-20	13 mm	20 mm
888-9516-20	16 mm	20 mm
888-9519-20	19 mm	20 mm



Screw for Dental Plugs


Titanium, LeiStar, T10, 3mm length

888-9501-00

Screwdriver Shaft for Dental Plugs

T30, AO Connection, Stainless Steel, 190mm in length


128-3019-00



Screwdriver Shaft for Dental Plug Screws

T10, AO Connection, Stainless Steel, 190mm in length

128-1019-00



Screwdriver T-Handle

Screwdriver Handle

Silicone, AO-Connection
sterilizable up to 134°C / 273°F



Product Code	Description
247-0102-00	standard, compatible with 2.7/3.5

Screwdriver T-Handle

Silicone, AO, not cannulated
sterilizable up to 134°C / 273°F



247-0101-00

